

Pasteur pipettes, Soda-lime glass Brand

www.huberlab.ch



Product Data

Made of Soda-lime glass

- ISO 7712
- Long drawn-out capillary tip
- Suction end with constriction for a cotton plug
- Volume: 2 ml
- Outer-Ø: 7 mm

Order unit:

1 unit with 250 pieces

additional packing information:

overpack 4 unit

| Item No. | length | execution | sterile | Price |
|-----------|--------|------------|---------|-----------|
| 1.1126.01 | 145 mm | short form | - | CHF 19.60 |
| 1.1127.01 | 225 mm | long form | - | CHF 22.21 |

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 15.05.2024 at 02:22 oclock.