


## Cold trap DURAN



Product in eShop

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Separation](#) / [Distillation & Extraction](#) / [Glass components & accessories](#) / [Coolers](#) / Cold trap DURAN

### Product Data

Made of borosilicate glass

Item No.	type no.	length	Outer dia.	external diameter	Price
5.6918.01	with half dip tube	220 mm	30 mm	olive 10 mm	CHF 144.05
5.6918.02	with completely dip tube	220 mm	30 mm	olive 10 mm	CHF 149.80
5.6918.09	2-tlg. mit seitlichen Oliven	---	---	Olive 10 mm mm	CHF 225.88

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 04.05.2024 at 16:29 oclock.

[www.huberlab.ch](http://www.huberlab.ch)